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- Guideline on Validation of the *Limulus* Amebocyte Lysate Test as an End-Product Endotoxin Test for Human and Animal Parenteral Drugs, Biological Products, and Medical Devices. U.S. Department of Health and Human Services, Public Health Service, Food and Drug Administration, December 1987.
- Bacterial Endotoxins Test, United States Pharmacopeia (current revision), United States Pharmacopoeial Convention, Rockville, MD.
- Bacterial Endotoxins, European Pharmacopoeia (current revision), European Pharmacopoeia Commission, Strasbourg, France.
- Bacterial Endotoxins Test, Japanese Pharmacopoeia (current revision), Tokyo, Japan.
- Lindsay, G. K., P. F. Roslansky, and T. Novitsky. 1989. Single-Step, Chromogenic Limulus Amebocyte Lysate Assay for Endotoxin J. Clinic. Microbiol. 27:947-951.
- Howell, W. H. 1885. Observations upon the chemical composition and coagulation of the blood of *Limulus polyphemus*, *Callinectes hastatus*, and *Cucumaria* sp. Johns Hopkins University Circular 43:4-5.
- Bang, F. B. 1953. The toxic effect of a marine bacterium on *Limulus* and the formation of blood clots. Biol. Bull. (Woods Hole, MA) 105:361-362.
- Levin, J., and F. B. Bang. 1964. A description of cellular coagulation in the *Limulus*. Bull. Johns Hopkins Hosp. 115:337-345.
- Levin, J., and F. B. Bang. 1968. Clottable protein in *Limulus*: its localization and kinetics of its coagulation by endotoxin. Thromb. Diath. Haemorrh. 19:186-197.
- Levin, J. and F. B. Bang. 1964. The role of endotoxin in the extracellular coagulation of *Limulus* blood. Bull. Johns Hopkins Hosp. 115:265-274.
- Nakamura, S., T. Morita, S. Iwanaga, M. Niwa, and K. Takahashi. 1977. A sensitive substrate for the clotting enzyme in horseshoe crab hemocytes. J. Biochem. 81:1567-1569.
- Novitsky, T. J. 1984. Discovery to Commercialization: The blood of the horseshoe crab. Oceanus 27:13-18.
- Progress in Clinical and Biological Research Vol. 231. Detection of Bacterial Endotoxins with the *Limulus* Amebocyte Lysate Test. 1987. Watson, S. W., J. Levin, and T. J. Novitsky (eds), Alan R. Liss, Inc., NY.
- Gould, M. C. Endotoxin in vertebrate cell culture: Its measurement and significance, p. 125- 136. In: Uses and Standardization of Vertebrate Cell Cultures, In Vitro Monograph number 5. 1984. Tissue Culture Association, Gaithersburg, MD.
- Ebner, C., D. Kraft, F. Prasch, R. Steiner, and H. Ebner. 1992. Type I allergy induced by Limulus Amoebocyte Lysate (LAL). Clinical and Experimental Allergy 22:417-419.
- Tsuji, K. and S. J. Harrison. 1978. Dry-heat destruction of lipopolysaccharide: Dry-heat destruction kinetics. Appl. Env. Microbiol. 36:710-714.
- Sweet, B. H. and J. F. Huxsoll. Depyrogenation by Dry Heat, ch. 12, p.101-108. In: Depyrogenation, Technical Report No. 7. 1985. Parenteral Drug Association, Inc., Philadelphia, PA.
- Transfusion and Infusion Assemblies and Similar Medical Devices, chapter <161>, USP, current revision, United States Pharmacopeial Convention, Rockville, MD.

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