

24-June-2013

RE: Pyros® EQS V. 1.2 SP1 Upgrade

Dear Valued Customer:

Associates of Cape Cod, Inc. (ACC) is pleased to announce the latest release of Pyros® EQS Endotoxin Quantitation Software (PEQS), version 1.2 SP1.

PEQS is software designed for use in *Limulus* Amebocyte Lysate (LAL) kinetic turbidimetric or chromogenic assays for the quantitation of gram-negative bacterial endotoxin. The application collects and processes data from the Pyros Kinetix® Flex and Pyros Kinetix® tube readers, stores the information in a database, and produces reports on test results and data trends. PEQS is designed to be compliant with 21 CFR Part 11.

As part of ACC's continuing efforts to incorporate customer feedback and refine our LAL data collection systems, we've made a number of changes in Pyros® EQS v1.2, including:

- Simplified and automated approach for installation and un-installation of the software
- Now supports Microsoft Windows 7 and Windows Server 2008 R2
- Now supports Oracle® Database 11g
- Improved stability and performance

If you wish to order these or any other products and services from ACC, please contact your local representative or email ACC at: custservice@acciusa.com. We appreciate your business, and will continue to provide you with leading products and services to meet your current and future needs.

Sincerely,



James P. Spitzer
Director of Marketing