



NORTHSTAR
BioProducts™

Associates of Cape Cod, Inc.

www.acciusa.com/bio

Chromatography


5.1

Specialty
Chromatography
Media

Endotoxin
Removal Resin

Affinity
Chromatography
Resins

Lectin Agarose Kits





Catalog No.

Product Description

Package Size

Specialty Chromatography Media

Endotoxin Removal Resins

Cellufine ET clean is poly(e-lysine) immobilized Cellufine (cellulose spherical beads). The beads bind and remove endotoxin from sample solution. The poly(e-lysine) is a microbial poly(amino acid) that consists of 30–35 lysine residues produced by *Streptomyces albulus*.

681-984-324	Cellufine™* ET clean L	10 mL
681-984-326	Cellufine™* ET clean L	50 mL
681-984-328	Cellufine™* ET clean L	500 mL
682-985-324	Cellufine™* ET clean S	10 mL
682-985-326	Cellufine™* ET clean S	50 mL
682-985-328	Cellufine™* ET clean S	500 mL

Affinity Chromatography Resins

800430-1	Chitooligo-Agarose	5 x 1 mL
800432-1	Fucose-Agarose	5 x 1 mL
800434-1	Lactose-Agarose	5 x 1 mL
800436-1	Maltose-Agarose	5 x 1 mL
800437-1	Melibiose-Agarose	5 x 1 mL
800441-1	Carbohydrate-Agarose Set	1 set

WGA, Lotus, PNA, CON-A, and GS-I are immobilized on agarose resin.

Lectin Agarose Kits

Our Lectin Agarose Kits are for Lectin affinity chromatography, purification of glycoproteins, structural analysis of sugar chain in glycoproteins.

800255-1	Lectin-Agarose Set III (<i>AAL, DSA, MAM & SSA</i>)	1 set
800257-1	Lectin-Agarose Set I (<i>Con A, LCA, RCA120 & WGA</i>)	1 set
800258-1	Lectin-Agarose Set II (<i>PHA-E4, PHA-L4, PNA & UEA-I</i>)	1 set

Please contact Customer Service for availability on bulk quantities.

Please refer to product insert for storage conditions.

**A trademark of Chisso Corporation*



NORTHSTAR
BioProducts™
Associates of Cape Cod, Inc.

Corporate Headquarters

124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536-4445
Tel: (888) 395-2221 or (508) 540-3444
Fax: (508) 540-8680
www.acciusa.com/bio

Customer Service: custservice@acciusa.com

Technical Service: techservice@acciusa.com

Contract Test Service: testservice@acciusa.com

Germany

Pyroquant Diagnostik GmbH
Opelstrasse 14
D-64546 Mörfelden-Walldorf, Germany
Tel: (49) 61 05-96 10 0
Fax: (49) 61 05-96 10 15
E-mail: service@pyroquant.de
www.pyroquant.de

United Kingdom

Associates of Cape Cod International, Inc.
Deacon Park, Moorgate Road
Knowsley
Liverpool L33 7RX
United Kingdom
Tel: (44) 151-547-7444
Fax: (44) 151-547-7400
E-mail: info@acciuk.co.uk
www.acciuk.co.uk