



124 Bernard E. Saint Jean Drive  
East Falmouth, MA 02536-4445

888.395.2221 / 508.540.3444  
custservice@acciusa.com

**Subject:** Important Updates for US Domestic Shipments re PyroSmart NextGen® Recombinant Cascade Reagent (rCR); Catalog Numbers PNG050-2 and PNG050-100

Date: November 21, 2025

Dear Valued Customer,

No response is required. You are receiving this notification as a customer of record for PyroSmart NextGen® Recombinant Cascade Reagent (rCR). We are writing to inform you of three important updates, effective January 2026, designed to enhance product clarity, improve accessibility, and support our sustainability initiatives.

**1. Ambient Temperature Shipping**

Beginning January 2026, PyroSmart NextGen® rCR will be shipped at ambient temperature. This change reduces cold packaging and plastic materials, lowering environmental impact while maintaining the same high standards of product quality and performance.

Note: There are no changes to end-user storage requirements; product should continue to be ***stored at 2–8°C upon arrival.***

**2. Product Labeling**

To align with current industry terminology and USP Chapter <86>, the product name on vial labels and packaging will be updated from “Recombinant LAL Reagent” to “Recombinant Cascade Reagent.” This is a naming clarification only and does not affect product composition or performance.

**3. Instructions for Use (IFU)**

Beginning January 2026, IFUs will no longer be provided in paper format. Each product will include a QR Code Card inside the packaging, linking directly to the product IFU web page. The reverse side of the card will include Customer Service contact information for users who cannot access the IFU electronically.

We appreciate your continued partnership and support of our sustainability efforts. If you have any questions or need further information, please contact us at [custservice@acciusa.com](mailto:custservice@acciusa.com) or your local representative. A digital version of this notification and other product updates is available at [www.acciusa.com](http://www.acciusa.com) under the “News & Notices” section.

Sincerely,

ACC – Product Support Team

CORP\_0811