

Lab Automation Within Reach

Automate the tedious. Focus on the critical.



Andrew
the pipetting robot



Single/Multichannel Pipetting

- Sartorius Bluetooth enabled, electronic pipettes
- Dynamic pipetting range of 0.2 µL to 10 mL
- Versatile Connected Tools to move plates and tubes across work deck

Flexible Working Deck

- Over 200 integrated consumables
- Perform complex actions with Connected Devices for heating/cooling, shaking, SPE, and more
- Switch between methods in less than 15 mins

Small Footprint

- Compact, benchtop design (min. L 1078 x W 414 x H 606.5 mm)
- Easy installation in a hood or biosafety cabinet with 2 full rows of consumables for a depth of 22"

OneLab Software

- Intuitive graphical drag-and-drop design interface
- OneLab protocol library with over 40 ready to use methods
- Compliant-ready, browser-based deployment options



Reproducibility
Automated pipetting to minimize error and increase efficiency.



Rapid Method Transfer
Develop simple methods in minutes, complex in days, and transfer methods in seconds, whether across the hall or the globe.



Promote Health
Reduce risk of repetitive strain injury and exposure to hazardous materials.

Automate your most tedious protocols with simple, ready-to-use actions in OneLab!

All your methods in one place. Accessible. Traceable.

OneLab
design & execute

Apply magnetic field	Flow-through extraction	Pipetting	Serial dilution	Concentration normalization	Shaking	Heat / cool
----------------------	-------------------------	-----------	-----------------	-----------------------------	---------	-------------

Easy to use, reproducible results, reduced errors....

- Boehringer Ingelheim Pharma, Select Science



Good results for everyday use. Allows full traceability of experiments... High-quality, user-friendly interface with modern design for great sampling and good return of investment.

- Andrew+ User, Select Science



Waters, Andrew Alliance, Andrew+, Pipette+, and OneLab are trademarks of Waters Corporation. All other trademarks are the property of their respective owners. ©2022 Waters Corporation. Produced in the U.S.A. December 2022 720007789EN KP-PDF

waters.com/andrewalliance

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™